



Cryo-EM Merit Badge Standard Operating Procedure

A Collaboration Among the NIH-Funded Cryo-EM National Centers

Common Screening with EPU Standard Operating Procedure

version 1.2

Authors: Yan Liu, Htet Khant

Approved Version Date: 08/2026

Created Date: 11/2022

1. Purpose

- 1.1. Screening cryo-EM grids with EPU software (ThermoFisher Scientific)

2. Definitions:

- 2.1. EPU is an automated data collection software from ThermoFisher Scientific (TFS).
- 2.2. TEMUI is the user interface of the microscope
- 2.3. TEM Image Analysis (TIA) must be running to ensure EPU is communicating with TEMUI
- 2.4. Imaging Presets in EPU:
 - 2.4.1 Atlas
 - a. Create a grid map with stitching image tiles at low dose
 - b. Optic and camera settings are relatively fixed for each microscope
 - 2.4.2 Grid Square
 - a. Depends on mesh size of the grid (can see one grid square)
 - b. Illumination covers the entire grid at a small dose rate.
 - c. Linear mode or counting mode (camera setting) with binning x2.
 - 2.4.3 Hole/Eccentric Height
 - a. Field of view covers 3~4 holes.
 - b. Illumination covers ~15µm at a small dose rate.
 - c. Linear mode or counting mode (camera setting) with binning x2.
 - 2.4.4 Data Acquisition
 - a. Beam diameter meets a parallel illumination and fringe free conditions.
 - b. Dose rate is within camera sensitivity range.
 - c. Beam edge doesn't reach adjacent targets.
 - d. Exposure time is not too long to achieve the required dose (~50 e/A², depends on sample).
 - e. Electron counting mode with dose fractionations.
 - 2.4.5 Auto Focus
 - a. Same as Data Acquisition with shorter exposure and 2x binning.
 - 2.4.6 Thon Ring
 - a. Same as Data Acquisition with 2 s exposure and 2x binning.
 - 2.4.7 Drift Measurement
 - a. Same as Data Acquisition except with shorter exposure and binning.
 - 2.4.8 Zero Loss
 - a. Same as Data Acquisition except with shorter exposure and binning.

3 Supplies & Equipment

- ☐ Microscope
- ☐ Software: TEMUI, TIA, EPU

4 Procedure:

4.4 Microscope and Software Checks

- 4.4.1 Microscope's states are within specifications such as vacuums, HT, emission and temperatures.
- 4.4.2 Stop Data Acquisition of the previous session on **EPU** and close column valves.
- 4.4.3 Cassette with grids have been loaded (see Loading SOP)



Cryo-EM Merit Badge Standard Operating Procedure

A Collaboration Among the NIH-Funded Cryo-EM National Centers


4.5 Making Grid Atlases with EPU

4.5.1 After grids are loaded, wait for all temperatures statuses to be green AND <-160C

Temperature Control		
— Status —		
All Nitrogen Temperature		
— Dewar levels —		
Autoloader	29 %	2 h 41 min
Column	54 %	7 h 07 min
— Temperatures —		
Docker	92.2 K	-181.0 °C
Holder	80.1 K	-193.1 °C
Cassette gripper	90.0 K	-183.2 °C
Cartridge gripper	91.8 K	-181.4 °C
Autoloader Dewar	78.2 K	-195.0 °C
Column Dewar	78.8 K	-194.4 °C

4.5.2 Select **Inventory** on the **Autoloader** panel and wait till the loaded grids have been detected correctly in the cassette

4.5.3 Open EPU software, click **Preview** to take an image using **Atlas** preset (Column valve should be open)




4.5.3.1 **IF** there is a beam shift, insert FluScreen. In **TEMUI Tune** tab, go to subsection **Direct Alignments** then select **Beam shift**. **DO NOT TOUCH GUN TILT AND GUN SHIFT**.
Proceed to 4.2.4 if no beam shift adjustment is needed



Cryo-EM Merit Badge Standard Operating Procedure


A Collaboration Among the NIH-Funded Cryo-EM National Centers



- 4.5.3.2 On the hand panel, use **MF-X** and **MF-Y** knobs to center beam over the green circle.
- 4.5.3.3 Click **Done** after adjustment. Press **R1** to lift the screen
- 4.5.3.4 Take another **Preview** to ensure the beam is centered
- 4.5.4 Select the **Atlas** tab in **EPU > Session Setup** task
- 4.5.5 Select **New Session**, and enter a Name for the session: Name should be yearmonthday-projectnumber-atlas, for example 20210414-CT21-atlas
Image format = MRC, Output folder = X, click **Apply**
- 4.5.6 In the **Screening** section, check the cartridges with grids and select **Start** to begin. It'll take 10-15 minutes to map each grid

4.6 Screening Grid Squares

- 4.6.1 After Atlasing, evaluate ice thickness and grid quality. A good grid should have at least three quarters of thin area




- 4.6.2 Load a grid of interest by clicking on the slot number and selecting **Load from sample**
- 4.6.3 When loading is done, pick a square (usually a thinner square). Right click on the square, then select **Move stage here**
- 4.6.4 Go back to **Preparation** tab, select **GridSquare** preset and take a **Preview**. Adjust the magnification so that you only see one square without much cutoff
- 4.6.5 Go to **AutoFunctions**, select **Hole/EucentricHeight** preset and run **Auto-eucentric by beam**




Cryo-EM Merit Badge Standard Operating Procedure

A Collaboration Among the NIH-Funded Cryo-EM National Centers

tilt to find the eucentric height



- 4.6.6 Go back to **Preparation** tab and take a preview at **Hole/EucentricHeight** preset
- 4.6.7 Move stage over a hole and capture at **DataAcquisition** preset (**DataAcquisition** preset allows you to set up optics for data collection); Make sure the pixel size is adjusted to your sample
- 4.6.8 **Autofocus (optional)**: this is performed before capturing an image at **DataAcquisition** preset. Move stage on carbon at **Hole/EucentricHeight**, and run **Autofocus** task with **Autofocus** preset under **AutoFunctions** tab
- 4.6.9 Evaluate ice thickness, particle density, homogeneity etc. Below is a good image:



- 4.6.10 Screen at least one hole on the same square and repeat 4.6.3-4.6.9 to screen other squares on the same grid. Load the next grid to screen if needed
- 4.6.10 Image Shift Calibration
 - 4.6.10.1 After screening and all presets are determined, Image Shift Calibration is required for each session
 - 4.6.10.2 On the atlas map, find a feature, such as a piece of ice contaminant on a square. Right click and select **Move stage here** to move the stage onto the square
 - 4.6.10.3 Set Presets as **GridSquare**, and **Preview**
 - 4.6.10.4 Move stage close to the feature and run **Auto-eucentric by beam tilt** with **Hole/EucentricHeight** preset under **Autofunctions**
 - 4.6.10.5 **Preview** at **Hole/EucentricHeight** and make sure you can find the feature and center it in




Cryo-EM Merit Badge Standard Operating Procedure

A Collaboration Among the NIH-Funded Cryo-EM National Centers

the image. If the feature isn't visible, you may need to insert the FluScreen and move the stage to find the feature

4.6.10.6 **Preview** at **DataAcquisition** preset and make sure the feature is visible. If not, you may need to insert the FluScreen and move the stage to find the feature

4.6.10.7 Select **Preparation** tab > **Calibrate Image Shifts** task, and click **Start Calibration**, and follow the instructions. Use **Proceed** to continue and double click on the feature to center it



5 Chemicals: N/A

6 Waste Disposal: N/A